

Introducing **NEW**

KERYflex™

An Innovative Solution to Nail Fungus and Other Conditions:

- Safe in-office application
- Restores the appearance of natural nail

Delivers Three Key Benefits to You:

- Models and restores the toenail
- Creates a smooth, natural-looking surface
- Is both flexible and durable



KERYflex™

Now Available only Through a KeryFlex-Trained Podiatrist

Do You Have Problems with Your Nails?

Restore their Appearance **Today!**

Before



After



Look as good as you feel!

Read more about this exciting, new nail procedure and ask your podiatrist about . . .

KeryFlex™ is a trademark of Pod-Advance, Inc., 21 Business Park Drive, Branford, CT 06405. KeryFlex is not a treatment for nail fungus and other nail dystrophies.

Do Your Toenails Look Like This?



For a variety of reasons, some people are much more likely to develop toenail fungus and other conditions that cause the nail to lose its natural appearance. While medicines are available to treat the fungus, these compounds may take some time to be effective and may never restore the nail to its normal appearance.

Your podiatrist has been recently trained in a new technique that improves the appearance of nails damaged by fungus and other conditions.

This new technique is the KeryFlex Nail Restoration System.

KeryFlex uses polymer resins and special activators to create a durable, yet flexible artificial nail when exposed to a certain frequency of ultraviolet light.

KeryFlex Improves the Appearance of Your Toenails



In as little as 15 minutes, your podiatrist can create a realistic, natural-looking nail in place of the damaged nail.

Your podiatrist does this, first, by removing most of the damaged nail, using a special nail clipper and file or dremel. After the old, damaged nail has been removed, the podiatrist will begin the process of restoring the appearance of your nail: he will apply a bonding agent to help the new nail attach to the remnant of your natural nail.

Next, the podiatrist will model KeryFlex Resin to your nail bed, cure the resin with a special frequency of ultraviolet light for two minutes, and then possibly reshape the nail further, if necessary.

KeryFlex Nails are Durable and Flexible

When the podiatrist has completed the procedure, you will leave the office with a natural-looking nail.

Women may put nail polish on this nail, and they can also apply nail polish remover to the nail – the nail will not be affected in the least.

Men and women can engage in their normal exercise routines after their toe nails have been remodeled by KeryFlex.